

For Research Use Only

TXNDC6 Polyclonal antibody

Catalog Number: 21235-1-AP



Basic Information

Catalog Number: 21235-1-AP	GenBank Accession Number: BC107054	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 260 ug/ml by Nanodrop;	GeneID (NCBI): 347736	Recommended Dilutions: WB: 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q86XW9	
Isotype: IgG	Full Name: thioredoxin domain containing 6	
Immunogen Catalog Number: AG15719	Calculated MW: 263 aa, 30 kDa	
	Observed MW: 30 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : K-562 cells,
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cell lysates were subjected to SDS PAGE followed by western blot with 21235-1-AP (TXNDC6 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.