

For Research Use Only

FOXJ3 Polyclonal antibody

Catalog Number: 21240-1-AP



Basic Information

Catalog Number: 21240-1-AP	GenBank Accession Number: BC151828	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 700 µg/ml by Nanodrop and 347 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 22887	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: forkhead box J3	
Isotype: IgG	Calculated MW: 622 aa, 69 kDa	
Immunogen Catalog Number: AG15732	Observed MW: 66-69 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : L02 cells, HepG2 cells
Species Specificity: human, mouse, rat	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



L02 cells were subjected to SDS PAGE followed by western blot with 21240-1-AP (FOXJ3 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.