

For Research Use Only

# SERPINB12 Polyclonal antibody

Catalog Number: 21289-1-AP



## Basic Information

<b>Catalog Number:</b> 21289-1-AP	<b>GenBank Accession Number:</b> BC103885	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 153 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 89777	
<b>Source:</b> Rabbit	<b>Full Name:</b> serpin peptidase inhibitor, clade B (ovalbumin), member 12	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 425 aa, 48 kDa	
<b>Immunogen Catalog Number:</b> AG15818		

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human, mouse

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data