

For Research Use Only

# ZNF480 Polyclonal antibody

Catalog Number: 21328-1-AP

1 Publications



## Basic Information

### Catalog Number:

21328-1-AP

### Size:

150ul , Concentration: 267 µg/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG15298

### GenBank Accession Number:

BC018783

### GeneID (NCBI):

147657

### UNIPROT ID:

Q8WV37

### Full Name:

zinc finger protein 480

### Calculated MW:

516 aa, 59 kDa

### Purification Method:

Antigen affinity purification

## Applications

### Tested Applications:

ELISA

### Cited Applications:

IHC

### Species Specificity:

human

### Cited Species:

human

## Notable Publications

Author	Pubmed ID	Journal	Application
Qingyun Zhu	37519705	Heliyon	IHC

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

