

For Research Use Only

SIRT4 Polyclonal antibody

Catalog Number: 21440-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

21440-1-AP

Size:

150ul , Concentration: 140 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16075

GenBank Accession Number:

BC109319

GeneID (NCBI):

23409

UNIPROT ID:

Q9Y6E7

Full Name:

sirtuin (silent mating type information regulation 2 homolog) 4 (S. cerevisiae)

Calculated MW:

314 aa, 35 kDa

Observed MW:

~41 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Min Liu	28871079	Nat Commun	WB
Shi-Ting Li	32694775	Nat Metab	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

