

For Research Use Only

Liprin Alpha 2 Polyclonal antibody, PBS Only

Catalog Number: 21442-1-PBS



Basic Information

Catalog Number:

21442-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12989

GenBank Accession Number:

BC104912

GeneID (NCBI):

8499

UNIPROT ID:

O75334

Full Name:

protein tyrosine phosphatase, receptor type, f polypeptide (PTPRF), interacting protein (liprin), alpha 2

Calculated MW:

1236 aa, 141 kDa

Observed MW:

150 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-PPFIA2 (IP:21442-1-AP, 4ug; Detection:21442-1-AP 1:500) with mouse brain tissue lysate 2120 ug. This data was developed using the same antibody clone with 21442-1-PBS in a different storage buffer formulation.