

For Research Use Only

IL-1F10 Polyclonal antibody

Catalog Number: 21530-1-AP



Basic Information

Catalog Number: 21530-1-AP	GenBank Accession Number: BC103966	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 267 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 84639	
Source: Rabbit	Full Name: interleukin 1 family, member 10 (theta)	
Isotype: IgG	Calculated MW: 152 aa, 17 kDa	
Immunogen Catalog Number: AG16154		

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Background Information

IL1F10 is a member of the interleukin 1 cytokine family. This cytokine is thought to participate in a network of interleukin 1 family members to regulate adapted and innate immune responses.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 21530-1-AP (IL1F10 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.