

For Research Use Only

MIRO1 Polyclonal antibody

Catalog Number: 21560-1-AP

1 Publications



Basic Information

Catalog Number:

21560-1-AP

Size:

150ul, Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15605

GenBank Accession Number:

BC060781

GeneID (NCBI):

55288

UNIPROT ID:

Q8IXI2

Full Name:

ras homolog gene family, member T1

Calculated MW:

618 aa, 71 kDa

Observed MW:

70-80 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Notable Publications

Author	Pubmed ID	Journal	Application
Hao Bai	38827170	Bioact Mater	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 21560-1-AP (MIRO1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.