

For Research Use Only

ZNF699 Polyclonal antibody, PBS Only

Catalog Number: 21562-1-PBS



Basic Information

Catalog Number:

21562-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15990

GenBank Accession Number:

BC109268

GeneID (NCBI):

374879

UNIPROT ID:

Q32M78

Full Name:

zinc finger protein 699

Calculated MW:

642 aa, 74 kDa

Observed MW:

74 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human brain tissue were subjected to SDS PAGE followed by western blot with 21562-1-AP (ZNF699 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21562-1-PBS in a different storage buffer formulation.