

For Research Use Only

Integrin Alpha 6 Polyclonal antibody

Catalog Number: 21709-1-AP



Basic Information

Catalog Number: 21709-1-AP	GenBank Accession Number: BC136456	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 3655	Recommended Dilutions: IHC 1:50-1:500
Source: Rabbit	Full Name: integrin, alpha 6	
Isotype: IgG	Calculated MW: 1073 aa, 119 kDa	
Immunogen Catalog Number: AG16335		

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human breast cancer tissue, human colon cancer tissue

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data