

For Research Use Only

# WNT10B Polyclonal antibody

Catalog Number: 21729-1-AP



## Basic Information

<b>Catalog Number:</b> 21729-1-AP	<b>GenBank Accession Number:</b> BC096353	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 233 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 7480	<b>Recommended Dilutions:</b> IHC 1:50-1:500 IF 1:10-1:100
<b>Source:</b> Rabbit	<b>Full Name:</b> wingless-type MMTV integration site family, member 10B	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 389 aa, 43 kDa	
<b>Immunogen Catalog Number:</b> AG16396		

## Applications

**Tested Applications:**

IF, IHC, ELISA

**Species Specificity:**

human, mouse, rat

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

**Positive Controls:**

IHC : human liver tissue,

IF : NIH/3T3 cells, HeLa cells

## Background Information

### Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1%BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data