

For Research Use Only

WNT10B Polyclonal antibody

Catalog Number: 21729-1-AP

1 Publications



Basic Information

Catalog Number:

21729-1-AP

Size:

150ul, Concentration: 233 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16396

GenBank Accession Number:

BC096353

GeneID (NCBI):

7480

UNIPROT ID:

O00744

Full Name:

wingless-type MMTV integration site family, member 10B

Calculated MW:

389 aa, 43 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:50-1:500

IF/ICC 1:10-1:100

Applications

Tested Applications:

IF, IHC, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

rabbit

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human liver tissue,

IF/ICC : NIH/3T3 cells, HeLa cells

Notable Publications

Author	Pubmed ID	Journal	Application
Zhiyuan Bao	38478424	Biol Reprod	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

