

For Research Use Only

ZNF76 Polyclonal antibody

Catalog Number: 21828-1-AP



Basic Information

Catalog Number: 21828-1-AP	GenBank Accession Number: BC000133	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 7629	Recommended Dilutions: WB 1:200-1:1000 IP 0.5-4.0 ug for IP and 1:200-1:1000 for WB
Source: Rabbit	Full Name: zinc finger protein 76 (expressed in testis)	
Isotype: IgG	Calculated MW: 62 kDa	
Immunogen Catalog Number: AG16267		

Applications

Tested Applications: IP, WB, ELISA	Positive Controls:
Species Specificity: human, mouse	WB: HeLa cells, IP: HeLa cells,

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data