

For Research Use Only

NRXN3 Polyclonal antibody

Catalog Number: 21849-1-AP



Basic Information

Catalog Number: 21849-1-AP	GenBank Accession Number: BC152457	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 420 µg/ml by Nanodrop;	GeneID (NCBI): 9369	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: neurexin 3	
Isotype: IgG	Calculated MW: 1061 aa, 117 kDa	
Immunogen Catalog Number: AG16293	Observed MW: 69 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : Neuro-2a cells,
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Neuro-2a cells were subjected to SDS PAGE followed by western blot with 21849-1-AP (NRXN3 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.