

For Research Use Only

ZNF541 Polyclonal antibody, PBS Only

Catalog Number: 21873-1-PBS



Basic Information

Catalog Number:

21873-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16192

GenBank Accession Number:

BC101052

GeneID (NCBI):

84215

UNIPROT ID:

Q9H0D2

Full Name:

zinc finger protein 541

Calculated MW:

1365 aa, 148 kDa

Observed MW:

148 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21873-1-AP (ZNF541 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 21873-1-PBS in a different storage buffer formulation.