

For Research Use Only

ZNF423 Polyclonal antibody

Catalog Number: 21878-1-AP



Basic Information

Catalog Number:

21878-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16350

GenBank Accession Number:

BC112315

GeneID (NCBI):

23090

UNIPROT ID:

Q2M1K9

Full Name:

zinc finger protein 423

Calculated MW:

1284 aa, 145 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : Raji cells, mouse brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Raji cells were subjected to SDS PAGE followed by western blot with 21878-1-AP (ZNF423 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.