

For Research Use Only

Semaphorin 3A Polyclonal antibody

Catalog Number: 21925-1-AP

3 Publications



Basic Information

Catalog Number:

21925-1-AP

Size:

150ul, Concentration: 247 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16477

GenBank Accession Number:

BC111416

GeneID (NCBI):

10371

UNIPROT ID:

Q14563

Full Name:

sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3A

Calculated MW:

771 aa, 89 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IF, IHC, ELISA

Cited Applications:

WB, IF

Species Specificity:

human

Cited Species:

human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Min Zhang	29370308	PLoS One	IF
Victor H Guaiquil	35427568	Exp Eye Res	IF
Ruyu Yang	38955028	Int Immunopharmacol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

