

For Research Use Only

ZNF560 Polyclonal antibody, PBS Only

Catalog Number: 21980-1-PBS



Basic Information

Catalog Number:

21980-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15921

GenBank Accession Number:

BC101040

GeneID (NCBI):

147741

UNIPROT ID:

Q96MR9

Full Name:

zinc finger protein 560

Calculated MW:

790 aa, 91 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

Zinc Finger Protein 560 (ZNF560) is a member of the Krüppel-type (C2H2) zinc finger protein family, characterized by multiple zinc finger domains that mediate DNA binding and transcriptional regulation.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HEK-293 cells using ZNF560 antibody (21980-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 21980-1-PBS in a different storage buffer formulation.