

For Research Use Only

IL-13 Polyclonal antibody

Catalog Number: 22085-1-AP

1 Publications



Basic Information

Catalog Number:

22085-1-AP

Size:

150ul , Concentration: 227 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17280

GenBank Accession Number:

BC096139

GeneID (NCBI):

3596

UNIPROT ID:

P35225

Full Name:

interleukin 13

Calculated MW:

146 aa, 16 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

rat

Notable Publications

Author	Pubmed ID	Journal	Application
Chuanfeng Fang	26722390	Int J Clin Exp Pathol	WB,IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

