

For Research Use Only

FAM167A Polyclonal antibody

Catalog Number: 22106-1-AP



Basic Information

Catalog Number: 22106-1-AP	GenBank Accession Number: BC104041	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 83648	Recommended Dilutions: WB 1:1000-1:6000
Source: Rabbit	UNIPROT ID: Q96KS9	
Isotype: IgG	Full Name: family with sequence similarity 167, member A	
Immunogen Catalog Number: AG17336	Calculated MW: 214 aa, 24 kDa	
	Observed MW: 25 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : U-87 MG cells,
Species Specificity: human	

Background Information

Family with sequence similarity 167, member A (FAM167A) has ubiquitously low expression in all tissues types throughout the body. In mouse it has a higher expression in the skin, B-cells, and spleen, but the same low expression in all other cell types. It's reported that FAM167A acts as an essential factor for BCR-ABL-independent TKI resistance in CML by activating the noncanonical NF-κB pathway.(PMID: 35241148)

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 22106-1-AP (FAM167A antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.