

For Research Use Only

SIDT2 Polyclonal antibody, PBS Only

Catalog Number: 22155-1-PBS



Basic Information

Catalog Number:

22155-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17775

GenBank Accession Number:

BC114522

GeneID (NCBI):

51092

UNIPROT ID:

Q8NBJ9

Full Name:

SID1 transmembrane family, member

2

Calculated MW:

832 aa, 94 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data