

For Research Use Only

IRE1/ERN1 Polyclonal antibody

Catalog Number: 22203-1-AP

5 Publications



Basic Information

Catalog Number:

22203-1-AP

Size:

150ul, Concentration: 560 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17504

GenBank Accession Number:

BC130405

GeneID (NCBI):

2081

UNIPROT ID:

O75460

Full Name:

endoplasmic reticulum to nucleus signaling 1

Calculated MW:

977 aa, 110 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human, mouse, bovine

Notable Publications

Author	Pubmed ID	Journal	Application
Chan-Hung Chen	39056589	Environ Toxicol	WB
Yexiao Yan	38821134	Free Radic Biol Med	WB
Jing Wu	37179046	Food Chem Toxicol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

