

For Research Use Only

MUC2 Polyclonal antibody

Catalog Number: 22307-1-AP

7 Publications



Basic Information

Catalog Number:

22307-1-AP

Size:

150ul, Concentration: 180 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002457

GeneID (NCBI):

4583

UNIPROT ID:

Q02817

Full Name:

mucin 2, oligomeric mucus/gel-forming

Calculated MW:

540 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human

Cited Species:

human, mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Ata Ur Rehman	36056976	AMB Express	IHC
Hidayat Ullah	39723032	Food Sci Nutr	WB,IHC,IF
Guanru Qian	39250508	PLoS Pathog	IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

