

For Research Use Only

VEGFR2 Polyclonal antibody

Catalog Number: 22457-1-AP

5 Publications



Basic Information

| | | |
|--|--|--|
| Catalog Number: 22457-1-AP | GenBank Accession Number: BC131822 | Purification Method: Antigen affinity purification |
| Size: 150ul , Concentration: 680 µg/ml by Bradford method using BSA as the standard; | GeneID (NCBI): 3791 | |
| Source: Rabbit | Full Name: kinase insert domain receptor (a type III receptor tyrosine kinase) | |
| Isotype: IgG | Calculated MW: 1356 aa, 152 kDa | |
| Immunogen Catalog Number: AG17934 | | |

Applications

Tested Applications:
ELISA

Cited Applications:
IF, WB

Species Specificity:
human

Cited Species:
human, rat, mouse

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|---------------|-----------|---------------------|-------------|
| Xian Wang | 36049428 | Phytomedicine | WB,IF |
| Aini Pan | 38016610 | Int J Biol Macromol | WB |
| Wei-Xiang Sun | 37480413 | Curr Med Sci | WB,IF |

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

