

For Research Use Only

# VEGFR2 Polyclonal antibody

Catalog Number: 22457-1-AP **5 Publications**



## Basic Information

<b>Catalog Number:</b> 22457-1-AP	<b>GenBank Accession Number:</b> BC131822	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul, Concentration: 680 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 3791	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P35968	
<b>Isotype:</b> IgG	<b>Full Name:</b> kinase insert domain receptor (a type III receptor tyrosine kinase)	
<b>Immunogen Catalog Number:</b> AG17934	<b>Calculated MW:</b> 1356 aa, 152 kDa	

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
WB, IF

**Species Specificity:**  
human

**Cited Species:**  
human, mouse, rat

## Notable Publications

Author	Pubmed ID	Journal	Application
Xian Wang	36049428	Phytomedicine	WB,IF
Aini Pan	38016610	Int J Biol Macromol	WB
Wei-Xiang Sun	37480413	Curr Med Sci	WB,IF

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data