

For Research Use Only

MPDZ/MUPP1 Polyclonal antibody

Catalog Number: 22583-1-AP



Basic Information

Catalog Number: 22583-1-AP	GenBank Accession Number: BC140793	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 450 ug/ml by Nanodrop;	GeneID (NCBI): 8777	Recommended Dilutions: WB: 1:500-1:2000
Source: Rabbit	UNIPROT ID: O75970	
Isotype: IgG	Full Name: multiple PDZ domain protein	
Immunogen Catalog Number: AG18222	Calculated MW: 2070 aa, 222 kDa	
	Observed MW: 250 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : SH-SY5Y cells, SK-N-SH cells, U-87 MG cells
Species Specificity: human	

Background Information

The PDZ domain is a protein motif that mediates protein-protein interactions. MUPP1(multi-PDZ domain protein-1) has 13 PDZ domains, each containing the structural characteristics of a prototypic PDZ motif, but lacks a catalytic domain. MUPP1 belongs to the PDZ domain protein family and is a structural homologant of PATJ (PALS1-associated tight-linking protein) (PMID: 23880463). MUPP1 is a scaffold protein and plays an important role in cytoskeletal organization, cell polarity, and cell proliferation(PMID: 25662616). Mutations in MUPP1 are associated with congenital hydrocephalus(PMID: 23240096).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 22583-1-AP (MUPP1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.