

For Research Use Only

CTAG1B Polyclonal antibody, PBS Only

Catalog Number: 22785-1-PBS



Basic Information

Catalog Number:

22785-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18784

GenBank Accession Number:

BC130362

GeneID (NCBI):

1485

UNIPROT ID:

P78358

Full Name:

cancer/testis antigen 1B

Calculated MW:

180 aa, 18 kDa

Observed MW:

20 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 22785-1-AP (CTAG1B antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 22785-1-PBS in a different storage buffer formulation.