

For Research Use Only

# Cytokeratin 28 Polyclonal antibody

Catalog Number: 22813-1-AP



## Basic Information

Catalog Number:

22813-1-AP

Size:

150ul, Concentration: 550 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18932

GenBank Accession Number:

BC148794

GeneID (NCBI):

162605

UNIPROT ID:

Q7Z3Y7

Full Name:

keratin 28

Calculated MW:

464 aa, 51 kDa

Observed MW:

51-53 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse skin tissue,

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 22813-1-AP (KRT28 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.