

For Research Use Only

ATP8B2 Polyclonal antibody

Catalog Number: 22913-1-AP



Basic Information

Catalog Number:

22913-1-AP

Size:

150ul , Concentration: 260 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19018

GenBank Accession Number:

BC030288

GeneID (NCBI):

57198

UNIPROT ID:

P98198

Full Name:

ATPase, class I, type 8B, member 2

Calculated MW:

1209 aa, 137 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

ELISA

Species Specificity:

human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

