

For Research Use Only

LAMC3-Specific Polyclonal antibody

Catalog Number: 22978-1-AP



Basic Information

| | | |
|--------------------------------------|--|---|
| Catalog Number: 22978-1-AP | GenBank Accession Number: NM_006059 | Purification Method: Antigen affinity purification |
| Size: 150ul | GeneID (NCBI): 10319 | |
| Source: Rabbit | Full Name: laminin, gamma 3 | |
| Isotype: IgG | Calculated MW: 1575 aa, 171 kDa | |
| Immunogen Catalog Number: AG19465 | | |

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

