

For Research Use Only

# STIL Polyclonal antibody

Catalog Number: 22991-1-AP

1 Publications



## Basic Information

Catalog Number:

22991-1-AP

Size:

150ul, Concentration: 480 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19124

GenBank Accession Number:

BC126223

GeneID (NCBI):

6491

UNIPROT ID:

Q15468

Full Name:

SCL/TAL1 interrupting locus

Calculated MW:

1287 aa, 143 kDa

Purification Method:

Antigen Affinity purified

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human

## Notable Publications

Author	Pubmed ID	Journal	Application
Meiling Wang	38823260	Transl Oncol	WB, IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

