

For Research Use Only

CENPV Polyclonal antibody, PBS Only

Catalog Number: 23005-1-PBS



Basic Information

Catalog Number:

23005-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19210

GenBank Accession Number:

BC137486

GeneID (NCBI):

201161

UNIPROT ID:

Q7Z7K6

Full Name:

centromere protein V

Calculated MW:

275 aa, 30 kDa

Observed MW:

30-33 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23005-1-AP (CENPV antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 23005-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using 23005-1-AP (CENPV antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 23005-1-PBS in a different storage buffer formulation.