

For Research Use Only

DSG4 Polyclonal antibody

Catalog Number: 23087-1-AP



Basic Information

Catalog Number:

23087-1-AP

Size:

150ul , Concentration: 1000 ug/ml by Nanodrop and 633 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19520

GenBank Accession Number:

BC132907

GeneID (NCBI):

147409

UNIPROT ID:

Q86SJ6

Full Name:

desmoglein 4

Calculated MW:

1040 aa, 114 kDa

Observed MW:

160 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse skin tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 23087-1-AP (DSG4 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.