

For Research Use Only

TMEM127 Polyclonal antibody, PBS Only

Catalog Number: 23142-1-PBS



Basic Information

Catalog Number:

23142-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19472

GenBank Accession Number:

BC039892

GeneID (NCBI):

55654

UNIPROT ID:

O75204

Full Name:

transmembrane protein 127

Calculated MW:

238 aa, 26 kDa

Observed MW:

26 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

IHC, IF/ICC, IP, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




IP result of anti-TMEM127 (IP:23142-1-AP, 3ug; Detection:23142-1-AP 1:500) with HeLa cells lysate 3000ug. This data was developed using the same antibody clone with 23142-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 23142-1-AP (TMEM127 antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 23142-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 23142-1-AP (TMEM127 antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 23142-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed NIH/3T3 cells using TMEM127 antibody (23142-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 23142-1-PBS in a different storage buffer formulation.