

For Research Use Only

FAM46B Polyclonal antibody

Catalog Number: 23149-1-AP

Featured Product

4 Publications



Basic Information

Catalog Number:

23149-1-AP

Size:

150ul, Concentration: 1000 ug/ml by Nanodrop and 647 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19523

GenBank Accession Number:

BC014160

GeneID (NCBI):

115572

UNIPROT ID:

Q96A09

Full Name:

family with sequence similarity 46, member B

Calculated MW:

424 aa, 47 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF, IHC

Species Specificity:

human

Cited Species:

human, rat

Notable Publications

Author	Pubmed ID	Journal	Application
Patrick A Fletcher	36093576	Glia	IHC
Tao Liang	30532005	Exp Mol Med	WB
Hongyang Sang	35116892	Transl Cancer Res	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

