

For Research Use Only

TMEM77 Polyclonal antibody, PBS Only

Catalog Number: 23182-1-PBS



Basic Information

Catalog Number:

23182-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19711

GenBank Accession Number:

BC047025

GeneID (NCBI):

128338

UNIPROT ID:

Q6UX65

Full Name:

transmembrane protein 77

Calculated MW:

266 aa, 30 kDa

Observed MW:

29 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells lysate was subjected to SDS PAGE followed by western blot with 23182-1-AP (TMEM77 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 23182-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer slide using 23182-1-AP (TMEM77 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 23182-1-PBS in a different storage buffer formulation.