

For Research Use Only

HDAC9 Polyclonal antibody, PBS Only

Catalog Number: 23328-1-PBS



Basic Information

Catalog Number:

23328-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19930

GenBank Accession Number:

BC152405

GeneID (NCBI):

9734

UNIPROT ID:

Q9UKV0

Full Name:

histone deacetylase 9

Calculated MW:

1011 aa, 111 kDa

Observed MW:

100 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

IHC, IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 23328-1-AP (HDAC9 Antibody) at dilution of 1:200 (under 40x lens). This data was developed using the same antibody clone with 23328-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed SH-SY5Y cells using HDAC9 antibody (23328-1-AP) at dilution of 1:400 and CoraLite@488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 23328-1-PBS in a different storage buffer formulation.