

For Research Use Only

TMTC1 Polyclonal antibody

Catalog Number: 23349-1-AP

1 Publications



Basic Information

Catalog Number:

23349-1-AP

Size:

150ul, Concentration: 1000 ug/ml by Nanodrop and 520 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19944

GenBank Accession Number:

BC042083

GeneID (NCBI):

83857

UNIPROT ID:

Q8IUR5

Full Name:

transmembrane and tetra tricopeptide repeat containing 1

Calculated MW:

882 aa, 99 kDa

Observed MW:

99 kDa, 65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

human, mouse

Cited Species:

human, mouse

Positive Controls:

WB : mouse bladder tissue, mouse colon tissue, mouse spleen tissue

Notable Publications

| Author | Pubmed ID | Journal | Application |
|--------------|-----------|-------------------|-------------|
| Na Yoon Paik | 39190102 | Cell Mol Life Sci | WB, IP |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse bladder tissue were subjected to SDS PAGE followed by western blot with 23349-1-AP (TMTC1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.