

For Research Use Only

TUFT1 Polyclonal antibody, PBS Only

Catalog Number: 23385-1-PBS

Featured Product



Basic Information

Catalog Number:

23385-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20020

GenBank Accession Number:

BC002933

GeneID (NCBI):

7286

UNIPROT ID:

Q9NNX1

Full Name:

tuftelin 1

Calculated MW:

44 kDa

Observed MW:

44 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 23385-1-AP (TUFT1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 23385-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 23385-1-AP (TUFT1 Antibody) at dilution of 1:200 (under 10x lens). This data was developed using the same antibody clone with 23385-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of HEK-293 cells using 23385-1-AP (TUFT1 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 23385-1-PBS in a different storage buffer formulation.