

For Research Use Only

ANXA2R Polyclonal antibody, PBS Only

Catalog Number: 23838-1-PBS



Basic Information

Catalog Number:

23838-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20807

GenBank Accession Number:

BC067873

GeneID (NCBI):

389289

UNIPROT ID:

Q3ZCQ2

Full Name:

chromosome 5 open reading frame 39

Calculated MW:

193 aa, 22 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 23838-1-AP (ANXA2R antibody at dilution of 1:50) (under 10x lens). This data was developed using the same antibody clone with 23838-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 23838-1-AP (ANXA2R antibody at dilution of 1:50) (under 40x lens). This data was developed using the same antibody clone with 23838-1-PBS in a different storage buffer formulation.