

For Research Use Only

NDUFC1 Polyclonal antibody, PBS Only

Catalog Number: 23842-1-PBS

Featured Product



Basic Information

Catalog Number:

23842-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20835

GenBank Accession Number:

BC107682

GeneID (NCBI):

4717

UNIPROT ID:

O43677

Full Name:

NADH dehydrogenase (ubiquinone) 1, subcomplex unknown, 1, 6kDa

Calculated MW:

76 aa, 9 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 23842-1-AP (NDUFC1 Antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 23842-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human heart tissue slide using 23842-1-AP (NDUFC1 Antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 23842-1-PBS in a different storage buffer formulation.