

For Research Use Only

RLN1 Polyclonal antibody

Catalog Number: 23852-1-AP



Basic Information

Catalog Number: 23852-1-AP	GenBank Accession Number: BC005956	Purification Method: Antigen Affinity purified
Size: 150ul , Concentration: 347 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 6013	Recommended Dilutions: IHC 1:50-1:500 IF 1:50-1:500
Source: Rabbit	Full Name: relaxin 1	
Isotype: IgG	Calculated MW: 185 aa, 21 kDa	
Immunogen Catalog Number: AG20880		

Applications

Tested Applications: IF, IHC, ELISA	Positive Controls:
Species Specificity: human	IHC : human prostate hyperplasia tissue, human testis tissue
Note-IHC: <i>suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0</i>	IF : HeLa cells,

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human prostate hyperplasia tissue slide using 23852-1-AP (RLN1 antibody) at dilution of 1:200 (under 10x lens).