

For Research Use Only

ATP5SL Polyclonal antibody

Catalog Number: 23866-1-AP



Basic Information

Catalog Number:

23866-1-AP

Size:

150ul, Concentration: 750 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20909

GenBank Accession Number:

BC013323

GeneID (NCBI):

55101

UNIPROT ID:

Q9NW81

Full Name:

ATP5S-like

Calculated MW:

158 aa, 19 kDa

Observed MW:

29-32 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : A549 cells, HeLa cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23866-1-AP (ATP5SL antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.