

For Research Use Only

OXNAD1 Polyclonal antibody

Catalog Number:23907-1-AP



Basic Information

Catalog Number:

23907-1-AP

Size:

150ul , Concentration: 500 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20998

GenBank Accession Number:

BC008322

GeneID (NCBI):

92106

UNIPROT ID:

Q96HP4

Full Name:

oxidoreductase NAD-binding domain containing 1

Calculated MW:

312 aa, 35 kDa

Observed MW:

35 kDa

Purification Method:

Antigen Affinity purified

Recommended Dilutions:

WB: 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, U-937 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23907-1-AP (OXNAD1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.