

For Research Use Only

# HEATR4 Polyclonal antibody

Catalog Number: 23944-1-AP



## Basic Information

<b>Catalog Number:</b> 23944-1-AP	<b>GenBank Accession Number:</b> BC047590	<b>Purification Method:</b> Antigen Affinity purified
<b>Size:</b> 150ul , Concentration: 600 µg/ml by Nanodrop and 293 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 399671	<b>Recommended Dilutions:</b> IHC 1:20-1:200
<b>Source:</b> Rabbit	<b>Full Name:</b> HEAT repeat containing 4	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 1026 aa, 117 kDa	
<b>Immunogen Catalog Number:</b> AG21056		

## Applications

<b>Tested Applications:</b> IHC, ELISA	<b>Positive Controls:</b> IHC : human testis tissue, human breast cancer tissue
<b>Species Specificity:</b> human	
<b>Note-IHC: suggested antigen retrieval with <i>TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0</i></b>	

## Background Information

### Storage


**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis slide using 23944-1-AP (HEATR4 Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis slide using 23944-1-AP (HEATR4 Antibody) at dilution of 1:50 (under 40x lens).