

For Research Use Only

# HEATR4 Polyclonal antibody, PBS Only

Catalog Number: 23944-1-PBS



## Basic Information

Catalog Number:

23944-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21056

GenBank Accession Number:

BC047590

GeneID (NCBI):

399671

UNIPROT ID:

Q86WZO

Full Name:

HEAT repeat containing 4

Calculated MW:

1026 aa, 117 kDa

Purification Method:

Antigen Affinity purified

## Applications

Tested Applications:

IHC, Indirect ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis slide using 23944-1-AP (HEATR4 Antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 23944-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human testis slide using 23944-1-AP (HEATR4 Antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 23944-1-PBS in a different storage buffer formulation.