

For Research Use Only

TLL2 Polyclonal antibody, PBS Only

Catalog Number: 23988-1-PBS



Basic Information

Catalog Number:

23988-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21143

GenBank Accession Number:

BC047411

GeneID (NCBI):

83887

UNIPROT ID:

Q9BWW7

Full Name:

tubulin tyrosine ligase-like family, member 2

Calculated MW:

592 aa, 67 kDa

Observed MW:

67 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

TLL2 belongs to the Tubulin Tyrosine Ligase-Like (TLL) family, binds to tubulin and catalyzes the polyglutamylation modification of tubulin, participating in the organization and dynamic regulation of the microtubule cytoskeleton.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 23988-1-AP (TLL2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 23988-1-PBS in a different storage buffer formulation.