

For Research Use Only

CNPY1 Polyclonal antibody

Catalog Number: 23990-1-AP



Basic Information

Catalog Number:

23990-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21160

GenBank Accession Number:

BC107598

GeneID (NCBI):

285888

UNIPROT ID:

Q3B7I2

Full Name:

canopy 1 homolog (zebrafish)

Calculated MW:

92 aa, 11 kDa

Purification Method:

Antigen Affinity purified

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

