

For Research Use Only

SAPK4 Polyclonal antibody

Catalog Number: 24016-1-AP



Basic Information

Catalog Number: 24016-1-AP	GenBank Accession Number: BC000433	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 187 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 5603	
Source: Rabbit	UNIPROT ID: O15264	
Isotype: IgG	Full Name: mitogen-activated protein kinase 13	
Immunogen Catalog Number: AG21240	Calculated MW: 38 kDa	
	Observed MW: 42 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data