

For Research Use Only

# ZMYM2 Polyclonal antibody, PBS Only

Catalog Number: 24112-1-PBS



## Basic Information

Catalog Number:

24112-1-PBS

Size:

100ug, Concentration: 1 mg/ml by  
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21370

GenBank Accession Number:

BC036372

GeneID (NCBI):

7750

UNIPROT ID:

Q9UBW7

Full Name:

zinc finger, MYM-type 2

Calculated MW:

1377 aa, 155 kDa

Observed MW:

155-180 kDa

Purification Method:

Antigen Affinity purified

## Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 24112-1-AP (ZMYM2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24112-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 24112-1-AP (ZMYM2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 24112-1-PBS in a different storage buffer formulation.



IP result of anti-ZMYM2 (IP:24112-1-AP, 4ug; Detection:24112-1-AP 1:600) with Raji cells lysate 1320 ug. This data was developed using the same antibody clone with 24112-1-PBS in a different storage buffer formulation.